

Cutimed® DebriClean

Fast and effective wound debridement

The innovative,
efficient and
easy to use
debridement
pad



Cutimed®

Choose Cutimed® DebrClean for fast and effective wound cleansing




Strong cleansing efficiency^{1,2}



Effective bacterial binding³



> 99 % biofilm removal in a laboratory² environment




The gently abrasive white monofilament loops are specially designed to remove and retain debris, biofilm and microorganisms



The abrasive blue looped fibers break up firm slough and remove it from the wound

Foam layer is integrated to absorb cleansing solution and offers a retention capacity for wound exudate



Ordering information Cutimed® DebrClean

REF no.	Size	Dressings per box
76170-00	9 x 13 cm	1 x 5

¹ JB4 40031810, In vitro evaluation of the cleansing effect of Cutimed® DebrClean and Cutisoft Cotton. August 2018, BSN medical data on file.

² JB4 40031908, In vitro evaluation of removal of viable biofilm by Cutimed® DebrClean and Cutisoft Cotton. August 2018, BSN medical data on file.

³ JB4 40031811, In vitro evaluation of the bacterial adhesion to Cutimed® DebrClean and Cutisoft Cotton. July 2018, BSN medical data on file.

BSN medical Australia & New Zealand
an Essity company

L6, 211 Wellington Road,
Mulgrave, Victoria Australia 3170
PO Box 337, Mount Waverley, Victoria 3149

Customer Service Australia
T 1300 276 633 F 1300 998 830
www.bsnmedical.com.au

Customer Service New Zealand
T 0508 276 111 F 0508 998 830
www.bsnmedical.co.nz



Cutimed®,
an Essity brand

ESS0293 10/2020